KIT FOTOMETRO QTEST - MANUALE D'USO UE

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LeadIX GmbH
Hans-Knöll-Straße 6
07745 Jena

QTEST PHOTOMETER KIT - EU USER MANUAL

Research Use Only/not for diagnostics. Assemble and operate according to these instructions.

1. INTENDED USE & SCOPE

The QTest Photometer is a portable, standalone device designed to read the color results of miraculix QTests. The EU version is supplied as a kit and must be assembled by the user before first use.

Important: The device does not include an internal battery, but it can be easily powered via any 5 V USB source (≥ 500 mA), such as a computer, laptop, or power bank.

2. YOUR CUSTOMER JOURNEY (START HERE)

- Assemble the device – **follow the online Assembly Guide (See section 10).**
- Connect the photometer to your computer via USB – **the device powers on and you can access the software.**
- Read the README – this is the quick-start instruction file included in the photometer software on your PC. It guides you through the essential steps, **so you can use your QTest Photometer correctly from the very first try.**
- Open the Step by Step Guide – **a user oriented walkthrough that explains each step and the reasoning behind it for consistent results and learning (See section 10).**
- Measure – **use the short instructions in this manual during daily operation.**

3. KIT CONTENTS


Your QTest Photometer Kit contains the following parts (see Assembly Guide for illustrations):

- Cuvette holder with light-blocking lid
- Adafruit TSL2591 Light Sensor Board (indicator LED disabled)
- IO Rodeo 595 nm LED board
- Adafruit PyBadge (pre-programmed)
- Enclosure, screws, nuts, spacers
- Mini speaker and battery holder
- USB A–micro-B cable
- Semi-micro cuvettes (≥ 800 μ L, disposable)

4. TECHNICAL SPECIFICATIONS

PARAMETER	SPECIFICATION
Optics	595 nm LED board (IO Rodeo Colorimeter)
Display	OLED with live output + menu navigation
Programs	LSD (blotter/drop), MDMA (Pills/Powders), Psilocybin
Cuvettes	Semi-micro (≥ 800 μ L)
Power supply	USB 5 V / ≥ 500 mA (no internal battery)

5. PHOTOMETER DISPLAY & CONTROLS

- On/Off Switch** → turns the device on/off
- Menu / Navigation buttons** ($\uparrow$ / $\leftarrow$) → navigate through modes
- Blank button (top right)** → calibrates with unused detection solution
- Display** → shows menu, absorbance values, and results
- Cuvette holder** → insert cuvette with triangle facing you
- MicroUSB port** → connect to USB power supply

6. OPERATING INSTRUCTIONS (QUICK USE)

Preparation

- Place the photometer on a flat, stable surface.
- Connect it to 5 V USB power supply (≥ 500 mA), such as a computer, laptop, or power bank.
- Allow the LED to warm up for 5–10 minutes.
- Always wear nitrile gloves, safety glasses, and a lab coat when handling reagents. Follow the safety instructions in the QTest manual carefully before performing the test.

Blanking (Baseline Calibration)

Before getting your QTest results you need to show the photometer what the “color noise” of your reagents looks like. The photometer is pre-calibrated for each QTest type. That means it already “knows” what zero looks like and automatically filters out the natural color of the reagents.

Still, we recommend performing your own blank – it’s the most reliable and scientifically correct way to ensure accurate results. Why? During storage or transport, reagents (especially the Psilocybin extraction solution) may slightly change color. Running a blank allows the photometer to recognize your current reagent color, so your measurements stay precise and trustworthy.

How to prepare a blank (recommended):

- Run a QTest in parallel, following the same steps as in the instructions, but without adding any sample – this becomes your blank.
- Pipette ≥ 800 μ L of the blank solution into a clean cuvette. Use clean tools only
- Insert the cuvette (triangle facing you) into the cuvette holder, close the black lid, and press Blank.
- Store the remaining blank in a clean, tightly sealed, labeled glass vial in the fridge (2–8 °C), protected from light.
 - Stable for up to 1 month.
 - ⚠ Do not return the used blank to the original reagent bottle. Dispose of it according to the instructions in the QTest manual.
 - 💡 For frequent testing, use Bulk Kits – extra reagent is included for blanks, so you don’t need to sacrifice a test.

One single QTest can provide enough solution to prepare 4–5 blanks when handled carefully. If you want to understand when blanking really makes a difference, visit our blog: “Why Blanking Matters” (QR code)

Measurement

- Perform the QTest according to the instructions provided with your kit.
- After the reaction is complete, pipette at least 800 μ L of the colored solution into a new, clean cuvette.
- Mix gently by pipetting up and down inside the cuvette to ensure homogeneous color.
- Insert the cuvette (triangle facing you) into the cuvette holder and place the lid.
- The photometer will display the quantitative result.

💡 Tip: You can measure several samples after one blank, as long as you remain within the same QTest type. If switching off or to another substance (e.g., from MDMA to LSD), a new blank is required.

7. STORAGE AND SHELF LIFE

- Store the photometer in a clean, dry environment, protected from dust and direct sunlight.
- Avoid extreme temperatures, humidity, or mechanical stress.
- Reagents are not included with this kit. Their shelf life and storage conditions are provided separately in the QTest manuals.
- Important:** Always remove the liquid from the cuvette immediately after measurement, as the solutions can damage the cuvette material.
- Never store any liquids inside the cuvettes or inside the photometer. Always clean and dry all parts thoroughly after use.

For information about storing the blank solution and its re-use stability, see section Blanking (Baseline Calibration).

8. SAFETY NOTES

- Always wear nitrile gloves, safety glasses, and a lab coat when handling reagents, following the safety requirements of the QTests (see QTest manual).
- Keep the device and reagents (when used) out of reach of children and animals.
- Avoid spills of reagents or acids on the photometer.
- Dispose of used cuvettes and reagents according to local regulations (see QTest manual).

9. DISCLAIMERS

- The QTest Photometer is provided for research purposes only.
- The device displays quantitative readouts for use with miraculix QTests.
- QTests are sold separately. Their use and limitations are described in their own manuals.
- Miraculix assumes no liability for misuse, incorrect assembly, or misinterpretation of results.
- Not CE-Certified

10. SUPPORT & RESOURCES

For correct use of your QTest Photometer, please consult the following resources:

- Assembly Guide (with photos): step-by-step online guide for building your device.
- README Instructions (inside software): quick start on blanking and measurement.
- Step-by-Step Guide (online): detailed scientific explanation and troubleshooting.

Assembly
Guide

README
Instructions

Step-by-Step
Guide